
Product Name: HNF 4 alpha Rabbit Monoclonal Antibody**Catalog #: AMRe85664**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Antigen Information


Gene Name	HNF 4 alpha
Alternative Names	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha
Gene ID	15378.0
SwissProt ID	P49698
Immunogen	Recombinant protein of mouse HNF-4-alpha

Background


Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and inTine.MiscellaneousBinds fatty acids.

Research Area

Image Data



Western blot analysis of HNF4alpha in mouse liver lysates using HNF 4 alpha antibody.



Western blot analysis of HNF4alpha in K562, C6, HeLa lysates using HNF4alpha antibody.