

Product Name: HPDL Rabbit Monoclonal Antibody
Catalog #: AMRe85673



Summary

Production Name	HPDL Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% BSA and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	HPDL
Alternative Names	4-hydroxyphenylpyruvate dioxygenase like; GLOXD1; 4-HPPD-L
Gene ID	84842.0
SwissProt ID	Q96IR7. A synthetic peptide of human HPDL

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 39 kDa

Background


May have dioxygenase activity.

Product Name: HPDL Rabbit Monoclonal Antibody
Catalog #: AMRe85673



Research Area

Image Data



Western blot analysis of HPDL in HeLa, Caco2, A549 lysates using HPDL antibody.

Note

For research use only.