

Product Name: HuC/HuD Rabbit Monoclonal Antibody**Catalog #: AMRe85684**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, IP 1:10-1:20
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Antigen Information


Gene Name	HuC/HuD
Alternative Names	Huc; mHuC; PLE21; 2600009P04Rik
Gene ID	15571/1996
SwissProt ID	Q60900/P26378
Immunogen	Recombinant protein of mouse HuC

Background

Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF mRNA. May also bind poly-A tracts via RRM 3. May be involved in neuronal differentiation and maintenance.

Research Area

Image Data



Western blot analysis of HuD/HuC in rat Brain lysates using HuC/HuD antibody.