

Product Name: Interferon gamma Receptor 1 Rabbit Monoclonal Antibody**Catalog #: AMRe85690**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 100 kDa

Antigen Information


Gene Name	Interferon gamma Receptor 1
Alternative Names	IFNGR1; Interferon gamma receptor 1; IFN-gamma receptor 1; IFN-gamma-R1; CDw119; CD antigen CD119
Gene ID	3459.0
SwissProt ID	P15260
Immunogen	A synthetic peptide of human IFNGR1

Background

Receptor for interferon gamma. Two receptors bind one interferon gamma dimer.

Research Area

Image Data



Western blot analysis of IFNGR1 in K562, C6, 3T3, HeLa lysates using Interferon gamma Receptor 1 antibody.