

Product Name: IGFBP2 Rabbit Monoclonal Antibody
Catalog #: AMRe85694



Summary

Production Name	IGFBP2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	IGFBP2
Alternative Names	IGFBP2; BP2; IBP2; Insulin-like growth factor-binding protein 2; IBP-2; IGF-binding protein 2; IGFBP-2
Gene ID	3485.0
SwissProt ID	P18065. Recombinant protein of human IGFBP2

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:10-1:20,IHC:1:50-1:100
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Background


Product Name: IGFBP2 Rabbit Monoclonal Antibody
Catalog #: AMRe85694



Inhibits IGF-mediated growth and developmental rates. IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.

Research Area

Image Data



Western blot analysis of IGFBP2 in T47D lysates using IGFBP2 antibody.

Note

For research use only.