
Product Name: IGFBP6 Rabbit Monoclonal Antibody**Catalog #: AMRe85695**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 35 kDa

Antigen Information


Gene Name	IGFBP6
Alternative Names	IBP6
Gene ID	3489.0
SwissProt ID	P24592
Immunogen	Recombinant protein of human IGFBP6

Background

IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.

Research Area

Image Data



Western blot analysis of IGFBP6 in K562, HeLa lysates using IGFBP6 antibody.