

Product Name: Integrin beta 5 Rabbit Monoclonal Antibody
Catalog #: AMRe85712

Summary

Production Name	Integrin beta 5 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Integrin beta 5
Alternative Names	Integrin beta-5
Gene ID	3693.0
SwissProt ID	P18084. Recombinant protein of human Integrin beta 5

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 88 kDa

Background


Product Name: Integrin beta 5 Rabbit Monoclonal Antibody
Catalog #: AMRe85712

Integrin alpha-V/beta-5 (ITGAV:ITGB5) is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand.
(Microbial infection) Integrin ITGAV:ITGB5 acts as a receptor for adenovirus type C.

Research Area

PI3K-Akt signaling pathway

Image Data



Western blot analysis of Integrin beta 5 in C6, CHO-K1, HeLa lysates using Integrin beta 5 antibody.

Note

For research use only.