
Product Name: LGALS3BP Rabbit Monoclonal Antibody**Catalog #: AMRe85752**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 65 kDa; Observed MW: 100 kDa

Antigen Information


Gene Name	LGALS3BP
Alternative Names	90K; M2BP; gp90; CyCAP; BTBD17B; MAC-2-BP; TANGO10B
Gene ID	3959.0
SwissProt ID	Q08380
Immunogen	Recombinant protein of human LGALS3BP

Background

Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

Research Area

Image Data



Western blot analysis of LGALS3BP in HT-1080, HepG2, A549 lysates using LGALS3BP antibody.