
Product Name: FABP1 Rabbit Monoclonal Antibody**Catalog #: AMRe85754**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
Molecular Weight	Calculated MW: 14 kDa; Observed MW: 14 kDa

Antigen Information


Gene Name	FABP1
Alternative Names	FABP 1; FABP1; FABP1; FABPL; FABPL_HUMAN; Fatty Acid Binding Protein 1; Fatty acid binding protein 1 liver; Fatty Acid Binding Protein; Fatty acid-binding protein 1; Fatty acid-binding protein; Fatty acid-binding protein liver; L FABP; L-FABP; liver; Liver-type fatty acid-binding protein.
Gene ID	2168.0
SwissProt ID	P07148
Immunogen	Recombinant protein of human liver FABP

Background

Binds free fatty acids and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

Research Area

Image Data



Western blot analysis of Liver FABP in mouse liver lysates using FABP1 antibody.