

Product Name: LMO2 Rabbit Monoclonal Antibody**Catalog #: AMRe85755**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, IP 1:10-1:20
Molecular Weight	Calculated MW: 18 kDa; Observed MW: 18 kDa

Antigen Information


Gene Name	LMO2
Alternative Names	TTG2; LMO-2; RBTN2; RHOM2; RBTN1
Gene ID	4005.0
SwissProt ID	P25791
Immunogen	A synthetic peptide of human LMO2

Background

Acts with TAL1/SCL to regulate red blood cell development. Also acts with LDB1 to maintain erythroid precursors in an immature state.

Research Area

Image Data



Western blot analysis of LMO2 in Jurkat, Hela, Raji lysates using LMO2 antibody.