

Product Name: Mesothelin Rabbit Monoclonal Antibody
Catalog #: AMRe85780



Summary

Production Name	Mesothelin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Mesothelin
Alternative Names	Pre-pro-megakaryocyte-potentiating factor; MPF
Gene ID	56047.0
SwissProt ID	Q61468. Recombinant protein of mouse Mesothelin

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 70 kDa; Observed MW: 40,70 kDa

Background


Membrane-anchored forms may play a role in cellular adhesion.

Product Name: Mesothelin Rabbit Monoclonal Antibody
Catalog #: AMRe85780



Research Area

Image Data



Western blot analysis of Mesothelin in 3T3, mouse lung lysates using Mesothelin antibody.

Note

For research use only.