
Product Name: MS2 Rabbit Monoclonal Antibody**Catalog #: AMRe85801**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 89 kDa; Observed MW: 89 kDa

Antigen Information


Gene Name	MS2
Alternative Names	MS2; CD156; CD156a
Gene ID	101.0
SwissProt ID	P78325
Immunogen	A synthetic peptide of MS2

Background


Possible involvement in extravasation of leukocytes.

Research Area

Image Data



Western blot analysis of MS2 in mouse pancreas lysates using MS2 antibody.



Western blot analysis of MS2 in K562 lysates using MS2 antibody.