
Product Name: MS4A14 Rabbit Monoclonal Antibody**Catalog #: AMRe85802**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 77 kDa; Observed MW: 77 kDa

Antigen Information


Gene Name	MS4A14
Alternative Names	Testis development protein NYD-SP21
Gene ID	84689.0
SwissProt ID	Q96JA4
Immunogen	Recombinant protein of human MS4A14

Background

Membrane-spanning 4-domains subfamily A member 14

Research Area

Image Data



Western blot analysis of MS4A14 in Jurkat, rat Brain, C6 lysates using MS4A14 antibody.