
Product Name: MTX3 Rabbit Monoclonal Antibody**Catalog #: AMRe85809**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Antigen Information


Gene Name	MTX3
Alternative Names	Metaxin-3
Gene ID	345778.0
SwissProt ID	Q5HYI7
Immunogen	A synthetic peptide of human MTX3

Background


Could function in transport of proteins into the mitochondrion.

Research Area

Image Data



Western blot analysis of MTX3 in rat heart lysates using MTX3 antibody.



Western blot analysis of MTX3 in mouse brain lysates using MTX3 antibody.