

Product Name: NDUFAF1 Rabbit Monoclonal Antibody**Catalog #: AMRe85832**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200,IP 1:10-1:20
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 35 kDa

Antigen Information


Gene Name	NDUFAF1
Alternative Names	CGI65; CIA30; Ndufaf1
Gene ID	51103.0
SwissProt ID	Q9Y375
Immunogen	A synthetic peptide of human NDUFAF1

Background


Chaperone protein involved in the assembly of the mitochondrial NADH:ubiquinone oxidoreductase complex (complex I).

Research Area

Image Data



Western blot analysis of NDUFAF1 in K562, HeLa lysates using NDUFAF1 antibody.



Immunocytochemistry analysis of NDUFAF1 (green) in LNcap using NDUFAF1 antibody, and DAPI (blue).

