
Product Name: NONO Rabbit Monoclonal Antibody**Catalog #: AMRe85862**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Antigen Information


Gene Name	NONO
Alternative Names	P54; NMT55; NRB54; P54NRB
Gene ID	4841.0
SwissProt ID	Q15233
Immunogen	A synthetic peptide of human NMT55/p54nrb

Background

DNA- and RNA binding protein, involved in several nuclear processes. Binds the conventional octamer sequence in double stranded DNA. Also binds single-stranded DNA and RNA at a site independent of the duplex site (By similarity). Involved in pre-mRNA splicing, probably as a heterodimer with SFPQ.

Research Area

Image Data



Western blot analysis of NMT55/p54nrb in 293T, MOLT4 lysates using NONO antibody.