
Product Name: NTPCR Rabbit Monoclonal Antibody**Catalog #: AMRe85875**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 21 kDa; Observed MW: 21 kDa

Antigen Information


Gene Name	NTPCR
Alternative Names	C1orf57; HCR-NTPase
Gene ID	84284.0
SwissProt ID	Q9BSD7
Immunogen	Recombinant protein of human Nucleoside triphosphate phosphohydrolase

Background

Has nucleotide phosphatase activity towards ATP, GTP, CTP, TTP and UTP. Hydrolyzes nucleoside diphosphates with lower efficiency.

Research Area

Image Data



Western blot analysis of Nucleoside triphosphate Phospho-hydrolase in K562, C6, 3T3, HeLa lysates using NTPCR antibody.