

Product Name: NUDT9 Rabbit Monoclonal Antibody**Catalog #: AMRe85877**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 39 kDa

Antigen Information


Gene Name	NUDT9
Alternative Names	NUDT10
Gene ID	53343.0
SwissProt ID	Q9BW91
Immunogen	A synthetic peptide of human NUDT9

Background


Hydrolyzes ADP-ribose (ADPR) to AMP and ribose 5'-phosphate.

Research Area

Image Data



Western blot analysis of NUDT9 in K562, HeLa lysates using NUDT9 antibody.



Immunocytochemistry analysis of NUDT9(green) in HeLa using NUDT9 antibody, and DAPI(blue)