
Product Name: Nuf2 Rabbit Monoclonal Antibody**Catalog #: AMRe85878**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

Antigen Information


Gene Name	Nuf2
Alternative Names	CDCA1; CT106; hNuf2; hNuf2R; hsNuf2; NUF2; NUF2R
Gene ID	83540.0
SwissProt ID	Q9BZD4
Immunogen	A synthetic peptide of human Nuf2

Background

Acts as a component of the essential kinetochore-associated NDC80 complex, which is required for chromosome segregation and spindle checkpoint activity. Required for kinetochore integrity and the organization of stable microtubule binding sites in the outer plate of the kinetochore.

Research Area

Image Data



Western blot analysis of Nuf2 in K562 lysates using Nuf2 antibody.