

Product Name: RAB9p40 Rabbit Monoclonal Antibody**Catalog #: AMRe85895**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Antigen Information


Gene Name	p40
Alternative Names	p40; RAB9p40; RABEPK; RGD1310612; RP11 65N13.1; RP23-446N16.2
Gene ID	10244.0
SwissProt ID	Q7Z6M1
Immunogen	Recombinant protein of human p40

Background

Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

Research Area

Image Data



Western blot analysis of p40 in K562, HeLa lysates using p40 antibody.