

Product Name: RRAS Rabbit Monoclonal Antibody
Catalog #: AMRe85903



Summary

Production Name	RRAS Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	RRAS
Alternative Names	p23
Gene ID	6237.0
SwissProt ID	P10301. A synthetic peptide of human RRAS

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Regulates the organization of the actin cytoskeleton (PubMed:16537651, PubMed:18270267). With OSPBL3, modulates

Product Name: RRAS Rabbit Monoclonal Antibody
Catalog #: AMRe85903




integrin beta-1 (ITGB1) activity (PubMed:18270267).

Research Area

TGF-beta signaling pathway, MAPK signaling pathway

Image Data



Western blot analysis of RRAS in HeLa lysates using RRAS antibody.

Note

For research use only.