
Product Name: FURIN Rabbit Monoclonal Antibody**Catalog #: AMRe86172**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	87 kD

Antigen Information


Gene Name	FURIN
Alternative Names	FUR; PACE; PCSK3; Furin; Dibasic-processing enzyme; Paired basic amino acid residue-cleaving enzyme; PACE;FURIN
Gene ID	5045.0
SwissProt ID	P09958
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human FURIN. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to FURIN

Research Area

Image Data



Western blot analysis of FURIN expression in HEK293T (A), HeLa (B), Jurkat (C), rat liver (D) whole cell lysates. (Predicted band size: 86 kD; Observed band size: 87 kD)