

Product Name: AKAP14 Rabbit Monoclonal Antibody**Catalog #: AMRe86173**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	22 kD

Antigen Information


Gene Name	AKAP14
Alternative Names	AKAP28; A-kinase anchor protein 14; AKAP-14; A-kinase anchor protein 28 kDa; AKAP 28; Protein kinase A-anchoring protein 14; PRKA14;AKAP14
Gene ID	158798.0
SwissProt ID	Q86UN6
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human AKAP14. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to AKAP14

Research Area

Image Data



Western blot analysis of AKAP14 expression in HEK293T (A), HepG2 (B), Jurkat (C), mouse brain (D), mouse liver (E), rat brain (F), rat liver (G) whole cell lysates. (Predicted band size: 22 kD; Observed band size: 22 kD)