
Product Name: XPO5 Rabbit Monoclonal Antibody**Catalog #: AMRe86174**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	136 kD

Antigen Information


Gene Name	XPO5
Alternative Names	KIAA1291; RANBP21; Exportin-5; Exp5; Ran-binding protein 21;Exportin 5
Gene ID	57510.0
SwissProt ID	Q9HAV4
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Exportin 5. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to Exportin 5

Research Area

Image Data



Western blot analysis of Exportin 5 expression in HEK293T (A), MCF7 (B), Jurkat (C), mouse brain (D), rat brain (E) whole cell lysates. (Predicted band size: 136 kD; Observed band size: 136 kD)