

Product Name: GABRA5 Rabbit Monoclonal Antibody**Catalog #: AMRe86175**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	70 kD

Antigen Information


Gene Name	GABRA5
Alternative Names	Gamma-aminobutyric acid receptor subunit alpha-5; GABA(A) receptor subunit alpha-5;GABRA5
Gene ID	2558.0
SwissProt ID	P31644
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human GABRA5. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to GABRA5

Research Area

Image Data



Western blot analysis of GABRA5 expression in 293T (A), A549 (B), K562 (C), mouse liver (D) whole cell lysates. (Predicted band size: 52 kD; Observed band size: 70 kD)