

Summary

Production Name	GTF2F2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	GTF2F2 RAP30; General transcription factor IIF subunit 2; ATP-dependent helicase GTF2F2;
Alternative Names	General transcription factor IIF 30 kDa subunit; Transcription initiation factor IIF subunit beta; TFIIF-beta; Transcription initiation factor RAP30; TFIIF RAP 30
Gene ID	2963.0
SwissProt ID	P13984. KLH-conjugated synthetic peptide encompassing a sequence within human TFIIF RAP 30. The exact sequence is proprietary.

Application

Dilution Ratio	WB:1:500-1:1000, ICC:1:50-1:200
Molecular Weight	28 kD

Product Name: GTF2F2 Rabbit Monoclonal Antibody
Catalog #: AMRe86177




Background

Recombinant rabbit monoclonal antibody to TFIIF RAP 30

Research Area

Image Data



Western blot analysis of TFIIF RAP 30 expression in A549 (A), HeLa (B), K562 (C) whole cell lysates. (Predicted band size: 28 kD; Observed band size: 28 kD)

Note

For research use only.