

Product Name: ARHGAP26 Rabbit Monoclonal Antibody**Catalog #: AMRe86178**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	100 kD

Antigen Information


Gene Name	ARHGAP26 GRAF; KIAA0621; OPHN1L; Rho GTPase-activating protein 26; GTPase regulator associated
Alternative Names	with focal adhesion kinase; Oligophrenin-1-like protein; Rho-type GTPase-activating protein 26;ARHGAP26
Gene ID	23092.0
SwissProt ID	Q9UNA1
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human ARHGAP26. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to ARHGAP26

Research Area

Image Data



Western blot analysis of ARHGAP26 expression in A549 (A), mouse muscle (B) whole cell lysates. (Predicted band size: 92 kD; Observed band size: 100 kD)