

Product Name: CLP1 Rabbit Monoclonal Antibody
Catalog #: AMRe86180



Summary

Production Name	CLP1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	CLP1
Alternative Names	HEAB; Polyribonucleotide 5'-hydroxyl-kinase Clp1; Polyadenylation factor Clp1; Polynucleotide kinase Clp1; Pre-mRNA cleavage complex II protein Clp1; CLP1
Gene ID	10978.0
SwissProt ID	Q92989. KLH-conjugated synthetic peptide encompassing a sequence within human CLP1. The exact sequence is proprietary.

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	48 kD

Product Name: CLP1 Rabbit Monoclonal Antibody
Catalog #: AMRe86180



Background

Recombinant rabbit monoclonal antibody to CLP1

Research Area

Image Data



Western blot analysis of CLP1 expression in HEK293T (A) whole cell lysates. (Predicted band size: 47; 40 kD; Observed band size: 48 kD)

Note

For research use only.