

Product Name: TNPO3 Rabbit Monoclonal Antibody**Catalog #: AMRe86188**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	100 kD

Antigen Information


Gene Name	TNPO3
Alternative Names	IPO12; Transportin-3; Importin-12; Imp12; Transportin-SR; TRN-SR;TNPO3
Gene ID	23534.0
SwissProt ID	Q9Y5L0
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human TNPO3. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to TNPO3

Research Area

Image Data



Western blot analysis of TNPO3 expression in HEK293T (A), Jurkat (B), mouse liver (C), mouse kidney (D), rat liver (E), rat kidney (F) whole cell lysates. (Predicted band size: 104 kD; Observed band size: 100 kD)