
Product Name: SRP19 Rabbit Monoclonal Antibody**Catalog #: AMRe86189**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	20 kD

Antigen Information


Gene Name	SRP19
Alternative Names	Signal recognition particle 19 kDa protein; SRP19;SRP19
Gene ID	6728.0
SwissProt ID	P09132
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human SRP19. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to SRP19

Research Area

Image Data



Western blot analysis of SRP19 expression in HEK293T (A), K562 (B), mouse liver (C), mouse kidney (D), mouse kidney (E), rat liver (F), rat kidney (G) whole cell lysates. (Predicted band size: 16 kD; Observed band size: 20 kD)