
Product Name: AGO3 Rabbit Monoclonal Antibody**Catalog #: AMRe86192**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	97 kD

Antigen Information


Gene Name	AGO3
Alternative Names	EIF2C3; Protein argonaute-3; Argonaute3; hAgo3; Argonaute RISC catalytic component 3; Eukaryotic translation initiation factor 2C 3; eIF-2C 3; eIF2C 3;eIF2C3
Gene ID	192669.0
SwissProt ID	Q9H9G7
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human eIF2C3. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to eIF2C3

Research Area

Image Data



Western blot analysis of eIF2C3 expression in HEK293T (A), mouse liver (B) whole cell lysates. (Predicted band size: 97; 71 kD; Observed band size: 97 kD)