
Product Name: COPS4 Rabbit Monoclonal Antibody**Catalog #: AMRe86193**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	50 kD

Antigen Information


Gene Name	COPS4
Alternative Names	CSN4; COP9 signalosome complex subunit 4; SGN4; Signalosome subunit 4; JAB1-containing signalosome subunit 4; COPS4
Gene ID	51138.0
SwissProt ID	Q9BT78
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human COPS4. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to COPS4

Research Area

Image Data



Western blot analysis of COPS4 expression in K562 (A) whole cell lysates. (Predicted band size: 46 kD; Observed band size: 50 kD)