

Product Name: MYH7B Rabbit Monoclonal Antibody
Catalog #: AMRe86197



Summary

Production Name	MYH7B Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	MYH7B
Alternative Names	KIAA1512; Myosin-7B; Antigen MLAA-21; Myosin cardiac muscle beta chain; Myosin heavy chain 7B cardiac muscle beta isoform; Slow A MYH14; MYH7B
Gene ID	57644.0
SwissProt ID	A7E2Y1.KLH-conjugated synthetic peptide encompassing a sequence within human MYH7B. The exact sequence is proprietary.

Application

Dilution Ratio	WB 1:500-1:1000, ICC 1:50-1:200
Molecular Weight	221 kD

Product Name: MYH7B Rabbit Monoclonal Antibody
Catalog #: AMRe86197




Background

Recombinant rabbit monoclonal antibody to MYH7B

Research Area

Image Data



Western blot analysis of MYH7B expression in mouse liver (A) whole cell lysates. (Predicted band size: 225; 20 kD; Observed band size: 221 kD)

Note

For research use only.