
Product Name: CSN1S1 Rabbit Monoclonal Antibody**Catalog #: AMRe86198**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	22 kD

Antigen Information


Gene Name	CSN1S1
Alternative Names	CASA; CSN1; Alpha-S1-casein;CSN1S1
Gene ID	1446.0
SwissProt ID	P47710
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human CSN1S1. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to CSN1S1

Research Area

Image Data



Western blot analysis of CSN1S1 expression in mouse liver (A), mouse kidney (B), rat liver (C) whole cell lysates. (Predicted band size: 21; 20 kD; Observed band size: 22 kD)