
Product Name: SRD5A2 Rabbit Monoclonal Antibody**Catalog #: AMRe86200**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	28 kD

Antigen Information


Gene Name	SRD5A2
Alternative Names	3-oxo-5-alpha-steroid 4-dehydrogenase 2; 5 alpha-SR2; SR type 2; Steroid 5-alpha-reductase 2; S5AR 2; Type II 5-alpha reductase;SRD5A2
Gene ID	6716.0
SwissProt ID	P31213
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human SRD5A2. The exact sequence is proprietary.

Background

Recombinant rabbit monoclonal antibody to SRD5A2

Research Area

Image Data



Western blot analysis of SRD5A2 expression in rat liver (A) whole cell lysates.
(Predicted band size: 28 kD; Observed band size: 28 kD)