
Product Name: SOX2 Rabbit Monoclonal Antibody**Catalog #: AMRe86215**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:100,IP 1:10-1:20
Molecular Weight	Calculated MW:34 kDa; Observed MW:34 kDa

Antigen Information

Gene Name	SOX2
Alternative Names	ANOP3; MCOPS3
Gene ID	6657
SwissProt ID	P48431
Immunogen	A synthetic peptide of human SOX2

Background

This intronless gene encodes a member of the SRY-related HMG-box (SOX) family of transcription factors involved in the regulation of embryonic development and in the determination of cell fate. The product of this gene is required for stem-cell

maintenance in the central nervous system, and also regulates gene expression in the stomach. Mutations in this gene have been associated with optic nerve hypoplasia and with syndromic microphthalmia, a severe form of structural eye malformation. This gene lies within an intron of another gene called SOX2 overlapping transcript (SOX2OT). [provided by RefSeq, Jul 2008]

Research Area

Image Data

NCCIT

kDa

180 -

130 -

100 -

70 -

55 -

40 -

35 - ←

25 -

Western blot analysis of extracts from NCCIT cells using SOX2 Rabbit Monoclonal Antibody at 1:1000.