

Product Name: HuD/HuC Rabbit Monoclonal Antibody
Catalog #: AMRe86259



Summary

Production Name	HuD/HuC Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, IHC-P, ICC/IF, FC, IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	HuD/HuC
Alternative Names	Huc; mHuC; PLE21; 2600009P04Rik
Gene ID	15571
SwissProt ID	Q60900.

Application

Dilution Ratio	WB: 1:1000 IHC-P: 1:500-1:2000 ICC/IF: 1:200-1:500 FC: 1:200-1:500 IP: 1:20-1:50
Molecular Weight	Calculated MW:40 kDa; Observed MW:40 kDa

Background


Product Name: HuD/HuC Rabbit Monoclonal Antibody
Catalog #: AMRe86259




Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF mRNA. May also bind poly-A tracts via RRM 3. May be involved in neuronal differentiation and maintenance.

Research Area

Image Data



Western blot detection of HuD/HuC in Rat Brain lysates using HuD/HuC antibody(1:1000 diluted).



Immunofluorescence analysis of Neuro-2a cells labelling HuD/HuC with AMRe86259.

Note

For research use only.