
Product Name: HSPA4 Rabbit Monoclonal Antibody**Catalog #: AMRe86297**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:100-1:200
Molecular Weight	Calculated MW:94 kDa; Observed MW:110 kDa


Antigen Information

Gene Name	HSPA4
Alternative Names	RY; APG-2; HSPH2; hsp70; hsp70RY; HEL-S-5a; HS24/P52
Gene ID	3308
SwissProt ID	P34932
Immunogen	A synthetic peptide of human HSPA4


Background

Research Area

Image Data



Western blot analysis of extracts from K562 cells using HSPA4 Rabbit Monoclonal Antibody at 1:1000.



Immunofluorescence analysis of HeLa cells labelling HSPA4 with AMRe86297.