

Product Name: FOXF1 Rabbit Monoclonal Antibody
Catalog #: AMRe86330



Summary

Production Name	FOXF1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	FOXF1
Alternative Names	FKHL5; ACDMPV; FREAC1
Gene ID	2294
SwissProt ID	Q12946.

Application

Dilution Ratio	WB: 1:1000-1:5000 ICC/IF: 1:200-1:500 FC: 1:100-1:500
Molecular Weight	Calculated MW:40 kDa; Observed MW:40 kDa

Background

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The


Product Name: FOXF1 Rabbit Monoclonal Antibody
Catalog #: AMRe86330



specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development. [provided by RefSeq, Jul 2008]

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using AMRe86330 at 1:1000.

Note

For research use only.