

Product Name: Apolipoprotein A II Rabbit Monoclonal Antibody
Catalog #: AMRe86331



Summary

Production Name	Apolipoprotein A II Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	Apolipoprotein A II
Alternative Names	apoAII; Apo-All; ApoA-II
Gene ID	336
SwissProt ID	P02652.

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	Calculated MW:11 kDa; Observed MW:11 kDa

Background


Product Name: Apolipoprotein A II Rabbit Monoclonal Antibody
Catalog #: AMRe86331



This gene encodes apolipoprotein (apo-) A-II, which is the second most abundant protein of the high density lipoprotein particles. The protein is found in plasma as a monomer, homodimer, or heterodimer with apolipoprotein D. Defects in this gene may result in apolipoprotein A-II deficiency or hypercholesterolemia. [provided by RefSeq, Jul 2008]

Research Area

Image Data



Western blot analysis of extracts from Human plasma tissue using AMRe86331 at 1:1000.

Note

For research use only.