

Product Name: BOB1 Rabbit Monoclonal Antibody**Catalog #: AMRe86387**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	Calculated MW:27 kDa; Observed MW:35 kDa

Antigen Information


Gene Name	BOB1
Alternative Names	BOB1; OBF1; OCAB; OBF-1
Gene ID	5450
SwissProt ID	Q16633
Immunogen	A synthetic peptide of human BOB1

Background

Transcriptional coactivator that specifically associates with either OCT1 or OCT2

Research Area

Image Data



Western blot analysis of extracts from Ramos cells using BOB1 Rabbit Monoclonal Antibody at 1:1000.