
Product Name: MLLT6 Rabbit Monoclonal Antibody**Catalog #: AMRe86483**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	Calculated MW:112 kDa; Observed MW:112 kDa

Antigen Information


Gene Name	MLLT6
Alternative Names	AF17
Gene ID	4302
SwissProt ID	P55198
Immunogen	Recombinant protein of human MLLT6

Background


ALL1-fused gene from chromosome 17 protein

Research Area

Image Data



Western blot detection of MLLT6 in HeLa, A549, HL-60, U2OS cell lysates using MLLT6 antibody (1:1000 diluted).



Western blot analysis of extracts from T24, C2C12 cells and Rat brain tissue using AMRe86483 at 1:1000.