

Product Name: JunB Rabbit Monoclonal Antibody**Catalog #: AMRe86519**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,ICC/IF 1:100-1:200,IP 1:10-1:100
Molecular Weight	Calculated MW:36 kDa; Observed MW:42,43 kDa

Antigen Information


Gene Name	JunB
Alternative Names	AP-1
Gene ID	3726
SwissProt ID	P17275
Immunogen	A synthetic peptide of human JunB

Background

Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.


Research Area

Image Data



Western blot detection of JunB in K562,C6,3T3,HeLa cell lysates using JunB antibody(1:1000 diluted).

C6



Western blot analysis of extracts from C6 cells using AMRe86519 at 1:1000.