

Product Name: Uteroglobin Rabbit Monoclonal Antibody**Catalog #: AMRe86538**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	Calculated MW:10 kDa; Observed MW:16 kDa

Antigen Information

Gene Name	Uteroglobin
Alternative Names	UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP
Gene ID	7356
SwissProt ID	P11684
Immunogen	A synthetic peptide of human Uteroglobin


Background

This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic

ligands. Defects in this gene are associated with a susceptibility to asthma. [provided by RefSeq, May 2010]

Research Area

Image Data



Western blot analysis of extracts from A549 cells using Uteroglobin Rabbit Monoclonal Antibody at 1:1000.