
Product Name: SAA4 Rabbit Monoclonal Antibody**Catalog #: AMRe86572**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:100,FC 1:20-1:50
Molecular Weight	Calculated MW:15 kDa; Observed MW:15 kDa

Antigen Information

Gene Name	SAA4
Alternative Names	CSAA; C-SAA
Gene ID	6291
SwissProt ID	P35542
Immunogen	A synthetic peptide of human SAA4

Background

Major acute phase reactant

Research Area

Image Data

Human plasma

kDa
180 -
130 -
100 -
70 -
55 -
40 -
35 -
25 -



Western blot analysis of extracts from Human plasma tissue using SAA4 Rabbit Monoclonal Antibody at 1:1000.