

Product Name: EGR2 Rabbit Monoclonal Antibody**Catalog #: AMRe86578**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:2000,FC 1:200-1:500
Molecular Weight	Calculated MW:50 kDa; Observed MW:50 kDa

Antigen Information


Gene Name	EGR2
Alternative Names	Egr-2; Zfp-6; Krox20; NGF1-B; Zfp-25; Krox-20
Gene ID	13654
SwissProt ID	P51774
Immunogen	A synthetic peptide of mouse EGR2

Background

Sequence-specific DNA-binding transcription factor (By similarity)

Research Area

Image Data



Western blot analysis of extracts from Mouse thymus tissue using EGR2 Rabbit Monoclonal Antibody at 1:1000.