
Product Name: Pbx3 Rabbit Monoclonal Antibody**Catalog #: AMRe86580**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	Calculated MW:47 kDa; Observed MW:47 kDa

Antigen Information


Gene Name	Pbx3
Alternative Names	Pre-B-cell leukemia transcription factor 3
Gene ID	5090
SwissProt ID	P40426
Immunogen	A synthetic peptide of human Pbx3

Background

Transcriptional activator that binds the sequence 5'-ATCAATCAA-3'.

Research Area

Image Data



Western blot analysis of extracts from Mouse brain tissue using Pbx3 Rabbit Monoclonal Antibody at 1:1000.