
Product Name: GRK3 Rabbit Monoclonal Antibody**Catalog #: AMRe86633**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:20-1:50
Molecular Weight	Calculated MW:80 kDa; Observed MW:80 kDa

Antigen Information

Gene Name	GRK3
Alternative Names	BARK2; ADRBK2
Gene ID	157
SwissProt ID	P35626
Immunogen	A synthetic peptide of human GRK3


Background

The beta-adrenergic receptor kinase specifically phosphorylates the agonist-occupied form of the beta-adrenergic and related G protein-coupled receptors. Overall, the beta adrenergic receptor kinase 2 has 85% amino acid similarity with beta adrenergic

receptor kinase 1, with the protein kinase catalytic domain having 95% similarity. These data suggest the existence of a family of receptor kinases which may serve broadly to regulate receptor function. [provided by RefSeq, Jul 2008]

Research Area

Image Data



Western blot analysis of extracts from Rat muscle tissue using GRK3 Rabbit Monoclonal Antibody at 1:3000.