

Product Name: Dab1 Rabbit Monoclonal Antibody**Catalog #: AMRe86673**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	Calculated MW:64 kDa; Observed MW:80 kDa

Antigen Information


Gene Name	Dab1
Alternative Names	scm; scr; yot; AI956902; scrambler; C630028C02Rik
Gene ID	13131
SwissProt ID	P97318
Immunogen	A synthetic peptide of human Dab1

Background

Adaptor molecule functioning in neural development. May regulate SIAH1 activity.

Research Area

Image Data



Western blot analysis of extracts from Mouse brain tissue using Dab1 Rabbit Monoclonal Antibody at 1:1000.