

Product Name: Pygopus 2 Rabbit Monoclonal Antibody**Catalog #: AMRe86707**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	Calculated MW:41 kDa; Observed MW:55 kDa

Antigen Information


Gene Name	Pygopus 2
Alternative Names	1190004M21Rik
Gene ID	90780
SwissProt ID	Q9BRQ0
Immunogen	A synthetic peptide of human Pygopus 2

Background


Involved in signal transduction through the Wnt pathway.

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using Pygopus 2 Rabbit Monoclonal Antibody at 1:1000.



Western blot analysis of extracts from A431 cells using AMRe86707 at 1:3000.